


Hemaray 51Vet

Auto Vet Hematology Analyzer



- WBC 5-part differentials, specialized for animal use;
- Minimal 15ul sample consumption;
- Whole Blood, Peripheral Blood, Pre-diluted blood;
- 29 parameters (include 4 research parameters);
- 3 reagents with Original QC and Calibrations;
- Integrated system with 10.4 inches color touch screen;
- Compact design





Specifications

• Parameterst

WBC, BAS#, BAS%, NEU#, NEU%, EOS#, EOS%,
LYM#, LYM%, MON#, MON%, RBC, HGB, HCT, MCV,
MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT,
P-LCR, P-LCC, ALY#, ALY%, LIC#, LIC%,
WBC Histogram, RBC Histogram, PLT Histogram,
2D DIFF Scatter gram*3,
3D DIFF Scatter gram*1

• Principles

Electrical impedance; Calorimetry
Laser scattering method

• Animals supporting

2 pre-defined animals(Cat, Dog)
4 user-defined animals

• Throughput

60 S/H

• Sample consumption

Whole Blood: 15ul (CBC+DIFF/ CBC)
Pre-diluted Blood: 20ul

• Data storage

100,000 results

Item	Linear range	Carry-over	CV(whole blood)
WBC	$(0\sim500) \times 10^9/L$	$\leq 0.5\%$	$< 2.0\%$
RBC	$(0\sim8.5) \times 10^{12}/L$	$\leq 0.5\%$	$< 1.5\%$
HGB	0-250g/L	$\leq 0.5\%$	$< 1.5\%$
HCT	0~67%	$\leq 0.5\%$	-
PLT	$(0\sim5000) \times 10^9/L$	$\leq 1\%$	$< 4.0\%$

Hemaray 51Vet

Auto Vet Hematology Analyzer



• Mode

CBC; CBC+DIFF

• QC

L-J, X-B, X-R

• Reagent

RD-50VDiluent; RL-5DV lyse
RL-50HV lyse
Concentrated Cleanser

• Interface

USB*4, LAN
Internal printer, internal card reader

• Working Environment

10°C~30°C; humidity $\leq 85\%$

• Power

100V-240V, 50Hz/60Hz, 400VA
Fuse: T3.15AL 250V, $\Phi 5 \times 20$

• Size

365mm*470mm*525mm
Weight: 35Kg



Rayto reserves the right to change specifications without prior notice.
P/N: ENG-Hemaray 51Vet 20190322

Rayto Life and Analytical Sciences Co.,Ltd.

Rayto Industrial Building, Shuangming Blvd South,
East Hi-Tech Park, Guangming New District,
518107 Shenzhen, P.R.China.
Tel:+86(755)88832350 Fax:+86(755)86168796
Email:info@rayto.com Http://www.rayto.com